BRAUN

ExactFit™ 5

Upper arm blood pressure monitor

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BP 6200

Product description 產品說明

BP 6200

This product is manufactured by Kaz Europe Sàrl under a license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany.

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
English

Intended use of Braun ExactFit 5

Braun upper arm blood pressure monitor has been developed for accurate and comfortable blood pressure measurements. The measuring accuracy of Braun's upper arm blood pressure monitor was tested at the time of manufacture and was proven by clinical research in accordance with ESH.

What you should know about blood pressure

Blood pressure constantly changes throughout the day. It rises sharply in the early morning and declines during the late morning. Blood pressure rises again in the afternoon and finally drops to a low level at night. Also, it may vary in a short period of time. Therefore, readings from successive measurements can fluctuate.

Blood pressure readings taken from a healthy 31-year-old male, measured at 5-minute intervals

Blood pressure measured in a doctor's office only provides a momentary value. Repeated measurements at home better reflect one's actual blood pressure values under everyday conditions.

Moreover, many people have a different blood pressure when they measure at home, because they tend to be more relaxed than when in the doctor's office. Regular blood pressure measurements taken at home can provide your doctor with valuable information on your normal blood pressure values under actual «everyday» conditions.

The World Health Organisation (WHO) has set up the following standard blood pressure values when measured at resting pulse.

Blood pressure (mmHg)	Normal values	Mild hypertension	Severe hypertension
SYS = systole (upper value)	up to 140	140-180	over 180
DIA = diastole (lower value)	up to 90	90-110	over 110

Intended use

- This monitor is intended for use in measuring the blood pressure and pulse rate of adults and uses the oscillometric method.
- This product is intended for household use only.
- This device is not intended as a substitute for regular check-ups by your doctor, please continue to visit your doctor on a regular basis for a professional reading.

Warning and precautions

- To ensure accurate measuring results, carefully read the complete use instructions.
- This product is intended for household use only. Keep product and batteries away from children.
- People suffering from cardiac arrhythmia, vascular constriction, arteriosclerosis in extremities, diabetes or users of cardiac pacemakers should consult their doctor before measuring their blood pressure themselves, since deviations in blood pressure values may occur in such cases.
- If you are under medical treatment or taking any medication, please consult your doctor first.
- The use of this blood pressure monitor is not intended as a substitute for consultation with your doctor.

Product description (See page 2-3, Fig. 1)

1. Start button
2. Memory button **M**
3. Date / time adjust button
4. Set button
5. Average button
6. User A / B switch
7. LCD display
8. Hose port
9. Connector
10. Arm cuff
11. Air hose
12. Battery compartment cover


Fig. 2

Inserting batteries (See Fig. 2-3)

- Remove the battery compartment cover at the bottom of the unit and insert 4 AA LR6 alkaline batteries with correct polarity (see symbol in the battery compartment).
- Note: always re-set date and time after replacing new batteries to make sure the measurement results are stored with correct date and time.

Only discard empty batteries. They should not be disposed of in the household waste, but at appropriate collection sites or at your retailer.


Fig. 3

Key rules for accurate blood pressure measurement

- Always take readings at the same time of day, ideally in the morning and evening, under the same conditions.
- Do not measure within 30 minutes after smoking or consuming coffee or tea.
- Take off wrist watch and jewelry before fitting the cuff on the measuring arm.
- While taking a measurement, sit down, relax, keep still, and do not move or speak.
- Wrap the cuff snugly around your arm. The cuff must be at heart level.
- Do not vibrate the unit during measurement, or the proper measurement will not be achieved.
- Perform measurement quietly in a relaxed position.
- Sit in a chair with your feet flat on the floor.
- Do not wrap the cuff over jacket or sweater sleeve or measurement cannot be done.
- Remove tight-fitting clothing from your left arm.

- Do not in any way twist the arm cuff.
- Do not inflate the monitor's cuff when it is not wrapped around the arm.
- Do not attempt to disassemble or change any parts of the monitor including the cuff.
- Do not drop the product or put it through strong impact.
- The device is not supposed to be used if your arm has any wound or injury.

Choosing the Right Cuff

For accurate measurement, it is important to choose the correct size cuff which best fits your upper arm. Choose the cuff size according to your arm circumference and make sure the bottom of the cuff is 2~3cm above your elbow.

- Small/Medium Cuff = 22~32cm arm circumference
- Large/XLarge Cuff = 32~42cm arm circumference

Applying the arm cuff

1. Insert the connector to the hose port of the device (Fig. 4).
2. Slide the end of arm cuff furthest from the tube through the metal ring to a loop. The smooth cloth should be on the inside of the cuff.


Fig. 4


Fig. 5


Fig. 6

3. If the cuff is located correctly, the velcro will be on the outside of the cuff and the metal ring will not touch the skin (Fig. 5).

4. Put your left arm through the cuff loop. The bottom of the cuff should be approximately (2~3 cm) above elbow. The tube should lie over the brachial artery on the inside of the arm (Fig. 6).


Fig. 7


Fig. 8



Fig. 9

5. Pull the cuff so that the top and bottom edges are tightened around your arm (Fig. 7).
6. When the cuff is positioned properly, press the Velcro firmly against the pile side of the cuff.
7. This cuff is suitable for use if the <<index>> mark falls within the <<ok range>> marked by two arrows when the cuff is tightened around your arm (Fig. 8).
8. Sit on a chair and place your arm on the table so that the cuff is at the same level as your heart (Fig. 9).

Display BP6200


Note: On BP6200, the backlight goes on when the unit is turned on and it will remain on until the device will be switched off.

Select mode

How to select user A /user B

- Make sure the Product is in power off mode.
- Slide the user switch to user A or user B, the current user mode will flash on the LCD display.


The WHO/ESH indicator to evaluate blood pressure data

This device has a blood pressure level indicator established according to WHO and the European Hypertension Society Guidelines (ESH) in 2007. For every measurement displayed on the screen, the cursor will indicate the blood pressure level with the corresponding color code, from green to red. You can use this classification daily to guide you to understand your blood pressure level. If you are really concerned by the classification level, you should consult your doctor.

Setting month, date and time

- Insert battery.
- The "year" 20__ will blink display, user can press adjust ▲ (3) to adjust year by increase step "1".
- User press set button ⌚ (4) to start the day setting, then the "day" will blink display, user can press adjust ▲ (3) to adjust day by increase step "1".
- User press set button ⌚ (4) to start the month setting, then the "month" will blink display, user can press adjust ▲ (3) to adjust month by increase step "1".
- User press set button ⌚ (4) to start the hour setting, then the "hour" will blink display, user can press adjust ▲ (3) to adjust hour by increase step "1".
- User press set button ⌚ (4) to start the minute setting, then the "minute" will blink display, user can press adjust ▲ (3) to adjust minute by increase step "1".
- User press set button ⌚ (4) to end all the date/time setting, all the blink are stopped.

Note: Holding the adjust button will scroll the value.

Taking a measurement

Wrap the cuff around the arm (see "applying the arm cuff" section above).

1. Sit upright on the chair to have a correct posture.
2. Press and release start button ⏻ (1), date/time and current user will be displayed.
3. Set user A/B switch (6) to A for user A or B for user B, LCD display will show user A or B symbol.
4. Press and release the start button ⏻ (1), all icons on display will be shown for 2 seconds. The device will adjust to zero automatically. The measuring blood pressure symbol will then flash on the display and the air pressure will automatically pump up to certain pressure level and start measurement.

Do not move or talk in the midst of taking blood pressure measurement.

5. After the air pressure has increased, the pulse is detected, the heart rate symbol ❤️ will start flashing.

6. LCD Display will show the results and WHO indicator arrow after measurement.

After taking blood pressure measurement, turn off the device by pressing the start button ⏻ (1) or automatically after 1 minute.

Memory function

Your blood pressure monitor can store 60 readings for both users each for the BP6200.


Storing measurement data

After each blood pressure measurement, the Systolic pressure, diastolic pressure, pulse rate and the time & date of specific day will be automatically stored. Memory #01 is always the most recent one. Once the memory is full, the oldest values will be overwritten.

Press memory button **M** (2) to review the stored data. The last memory data (sys/dia/pul) with measurement date/time, and WHO indicator will show on LCD. Press memory button **M** (2) again to show the previous data. Be aware that the correct user A or B is chosen.

Average function


Press the average button **WA** (5) to show the full day average of past 7 days on LCD. Press the average button second to show the morning average of past 7 days on LCD. If the result is with morning hypertension, the morning hypertension icon ☀️ will be shown. Press the average button third time to show the evening average of past 7 days on LCD. Press the average button fourth time to show the full day average of past 7 days again.

Erasing data

Make sure the Product is in power off mode.

Press and hold memory button **M** (2) for more than 5 seconds, LCD will blink display "dEL ALL" (if the "slide switch is in User A side, the user A icon will show "☺️" or "dEL ALL" (if slide switch is in User B side, the user B icon will show)).

Press memory button **M** (2) again, LCD will display "----" to mean all the stored data of corresponding user have been deleted.

Irregular heart beat detector

The appearance of this symbol  signifies that a certain pulse irregularity was detected during the measurement. Talking, moving, shaking or an irregular pulse during the measurement can result in the appearance of this icon. Usually this is not a cause for concern, however if the symbol appears often, we recommend you seek medical advice. The device does not replace a cardiac examination, but serve to detect pulse irregularities at an early stage.

Low battery indicator

When the low battery indicator  flashes on the display, it means the battery is low and the four batteries need to be replaced with alkaline LR6 (AA) batteries.

** After replacing batteries, BPM will go into time-setting mode automatically and show the last measurement time on screen. Please set the current date / time before taking the next measurement to get correct average result.

Storage and cleaning

- Always keep the unit in the carrying case after its use.
- Do not put the item directly under the sunlight, in high temperature or humid and dusty places.
- Do not store in extremely low (less than -20°C) or high (more than 60°C) temperature.
- Use a piece of cloth with water or mild cleansing agent to clean the case and then use a piece of dry cloth to wipe it dry. Use a piece of dry cloth to wipe the cuff when it is dirty.
- Do not use any strong cleansers to clean it.
- When the unit is not to be used for a long time, remove the batteries. (Batteries may leak or cause harm).
- Do not modify the device. NEVER open the device! This will make the manufacturer's warranty invalid.

Calibration

This device was designed and manufactured for a long service life, however it is generally recommended to have the device inspected once a year to ensure correct function and accuracy.

Please contact the authorized service centre located in your country. Note: The calibration checking is not a free service. Please contact the Authorized Service Centre to get a quotation before you send out the product.

This device is not intended as a substitute for regular check-ups by your doctor, please continue to visit your doctor on a regular basis for a professional reading.

Production date is given by the LOT located in the back of the device. The first 3 numeric digits after the LOT No represents the day of the year of manufacture. The next 2 numeric digits represent the last two numbers of the calendar year of manufacture and the letters at the end designate the manufacturer of the product. E.g. LOT 15614VTN this product was made on the day 156, year 2014 at manufacturer identifier VTN.

What to do if

Problem	Reason	Solution
Heart rate symbol 	Appears in the measurement condition and flashes when pulse is detected.	• Measurement in progress, remain quiet.
Low battery indicator 	Appears when the battery voltage is excessively low or the positions of the batteries are incorrect.	• Replace all four batteries with new ones. Insert the batteries in the correct positions. Be aware of +/- positions.
Measuring error 	Appears when the accurate blood pressure and pulse could not be obtained.	<ul style="list-style-type: none"> • Press "start/stop" button again and remeasure. • Check if cuff is wrapped according to instructions. • Check that there are no kinks in tube. • Check palm if exerting effort. • Check if talking or moving during measurement. • Check if posture is correct.
E1 shows	The cuff is not secure	• Refasten the cuff and then measure again.
E2 shows	The cuff is very tight	• Refasten the cuff and then measure again.
E10 or E11 shows	The monitor detected motion, talking or the pulse is too poor while measuring.	• Relax for a moment and then measure again.
E20 shows	The measurement process does not detect the pulse signal.	• Loosen the clothing on the arm and then measure again.
E21 shows	Measure incorrectly	• Relax for a moment and then measure again.
EE3 - EE15 shows	Error while measuring	• Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty in page13 for contact information and return instructions.

Specifications

Method of measurement	Oscillometric
Model number	BP 6200
Range of measurement	Pressure 0~300 mmHg Pulse 40~199 beats/minute
Accuracy	Pressure +/- 3 mmHg Pulse +/- 5 % Max.
Inflation	Deluxe automatic
Display	Liquid crystal display – systolic, diastolic, pulse rate Backlight Display
Sets of memory	60 sets per user
Cuff size	Small cuff = 22-32cm arm circumference Large cuff = 32-42cm arm circumference
Operating temperature	+10 °C ~ + 40 °C, less than 85 % R.H.
Storage temperature	-20 °C ~ +60 °C, less than 85 % R.H.
Operating atmospheric pressure	860-1060 hPa
Unit weight	Approximately 500gr (without batteries)
Power supply	Alkaline battery: 4 x AA (LR6) 1.5V
Battery life	300 times measurement
Auto power off	Whenever not used for 1 minute
Accessories	4 batteries, 2 arm cuffs with tube, instruction manual, pouch, travel bag
Service life	5 years

IMPORTANT

Read the operating instructions.

If device is not used within specified temperature, humidity and atmospheric pressure ranges the technical accuracy of the measurement cannot not be guaranteed.

Classification:

- Internally powered equipment
- Type BF equipment
- IP22: Protected against solid foreign objects of 12.5 mm diameter and greater.
Protected against vertically falling water drops when the device is tilted up to 15°
- Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide
- Continuous operation with short-time loading

Operating temperature

Storage temperature

Storage humidity

This product conforms to the provisions of the EC directive 93/42/EEC (Medical Device Directive). This device conforms to the following standards:

- EN 60601-1: 2006 + AC:2010: - General requirements for basic safety and essential performance
- EN 60601-1-2:2007 - Electromagnetic compatibility requirements and tests
- EN 60601-1-11:2010 - Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- EN 1060-1:1995 + A2:2009 – Non-invasive sphygmomanometers - general requirements
- EN 1060-3:1997 + A2:2009 - Non-invasive sphygmomanometers - Supplementary requirements for electro-mechanical blood pressure measuring systems.
- EN 1060-4:2004 - Non-invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.
- IEC60601-1-11:2010-General requirements for basic safety and essential performance
Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

Medical electrical equipment needs special precautions regarding EMC.
For detailed description of EMC requirements please contact an authorized local service centre

Portable and mobile RF communications equipment can affect medical electrical equipment.

Please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at your local retailer or at appropriate collection points provided in your country.

Warranty

This product comes with limited warranty commencing on the date of purchase. Within the warranty period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, by replacing the complete appliance.

This warranty is applicable only for the appliance supplied by the appointed distributor.

This warranty does not cover: damage due to improper use, normal wear or use as well as defects that have a negligible effect on the value or operation of the appliance. The warranty becomes void if repairs are undertaken by unauthorized persons and if original Braun parts are not used. To obtain service within the warranty period, please contact us.

HK Listing No.: 130467

Guidance and manufacturer's declaration – electromagnetic emissions		
The ME equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the ME equipment should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The ME equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class BF	Complies
Harmonic emissions IEC 61000-3-2	Not Applicable	The ME equipment is solely battery powered.
Voltage fluctuations / flicker emissions	Not Applicable	

Non-Life Support Equipment Separation Distance Calculation (3Vrms / 3V/m compliance)			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz in ISM bands $d = \left[\frac{3.5}{V_t} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_t} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_t} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

Guidance and manufacturer's declaration – electromagnetic immunity			
The ME is intended for use in the electromagnetic environment specified below. The customer or the user of the ME should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV Contact ±8kV Air	Complies	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2.5GHz	Complies	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3 V/m.
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	Not Applicable (no electrical cabling)	Interference may occur in the vicinity of equipment marked with the following symbol:  Separation distance calculation provided above. If a known transmitter is present the specific distance can be calculated using the equations.
Electrical fast transient IEC 61000-4-4	±2kV power line ±1kV I/O lines	Not Applicable	The ME equipment is solely battery powered.
Surge IEC 61000-4-5	±1kV differential ±2kV common	Not Applicable	
Power frequency magnetic field IEC 61000-4-8	3 A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interrupts and voltage variations on power supply input lines IEC 61000-4-11	>95% dip 0.5 cycle 60% dip 5 cycles 70% dip 25 cycles 95% dip 5 sec.	Not Applicable	The ME equipment is solely battery powered.

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- EN 60601-1:2006 + AC:2010: - General requirements for basic safety and essential performance
- EN 60601-1-2:2007 - Electromagnetic compatibility requirements and tests
- EN 60601-1-11:2010 - Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- EN 1060-1:1995 + A2:2009 - Non-invasive sphygmomanometers - general requirements
- EN 1060-3:1997 + A2:2009 - Non-invasive sphygmomanometers - Supplementary requirements for electro-mechanical blood pressure measuring systems.
- EN 1060-4:2004 - Non-invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.
- IEC60601-1-11:2010-General requirements for basic safety and essential performance-Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

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© 2014, Kaz Europe Sàrl